

StellARray™ Gene Expression System



Contents

1.0	Introduction	1
2.0	Principles	2
3.0	Procedure Overview	3
4.0	StellARray™ qPCR Protocol at a Glance	4
5.0	Plate Format	5
6.0	Materials Needed (not provided)	5
7.0	Recommendations for RNA Sample Preparation and Quality Control	6
8.0	cDNA Synthesis	7
9.0	Protocol for qPCR	7
10.0	Interpretation of Results	10
11.0	References	11
12.0	FAQs	11
13.0	Ordering Information	13
14.0	Related Products	13

1.0 Introduction

Gene expression profiling is widely regarded to be integral for advancing life science research and clinical diagnosis in the 21st century. It is increasingly recognized that subsets of genes (hundreds rather than thousands) provide the most valuable gene profile definitions of any specific biological process. The StellARray™ Gene Expression System is a quantitative polymerase chain reaction (qPCR)-based method and provides reliable profiling of biologically focused gene sets. The StellARray™ Gene Expression System is comprised of three parts:

1. GeneSieve™ Gene Query
2. StellARray™ qPCR Arrays
3. Global Pattern Recognition™ (GPR) Data Analysis Tool

The GeneSieve™ Query merges databases containing gene-centric biological information in a new and innovative way. GeneSieve™ Query can filter or “sieve” thousands of genes for a custom population of real-time PCR arrays that have a unifying biological theme, such as a particular disease or pathway of interest.

Each StellARray™ qPCR Array is compiled to provide the highest information content with a scientifically relevant context. Each qPCR array is a 96- or 384-well formatted, wet lab validated group of gene primer sets, focused on specific biological processes. StellARray™ qPCR Arrays facilitate accurate gene expression analysis, as well as analysis of copy number variations (CNV). Therefore, it is possible to reveal a two-dimensional barcode of samples, such as tumors relating gene copy number with transcription levels.

StellARray™ qPCR Arrays can be used with most common qPCR thermal cyclers (for compatibility with your thermal cycler, please check table in “Plate Formats”, page 5) and the SYBR® Green-based Reagents of your preference.

The Global Pattern Recognition™ (GPR) Analysis Tool allows gene expression analysis to be simple, fast and reliable. Data is simply uploaded to the GPR program and the software uniquely positions the expression level of each gene with respect to all genes within an experiment. This negates the assumption that a gene (normalizer) has an invariant expression level. Global Pattern Recognition™ Analysis Tool is unbiased in that it enables the experimental

data to define the invariant normalizer genes, not the experimenter. The use of any gene as a potential normalizer also maximizes the use of the limited real estate on a *StellARray*™ qPCR Array by eliminating the loss of wells used to contain housekeeping genes as potentially erroneously predefined normalizers.

2.0 Principles

The traditional approach to measure gene expression changes from real-time PCR data has been to normalize the results of a gene of interest with respect to one or a selection of housekeeping genes (e.g. reference or normalizer genes). The general assumption is that the level of expression of the normalizer gene(s) does not change in the context of the experiment and can be used to normalize the variability in RNA quantity between individual samples. By normalizing to a housekeeping gene, in theory, a magnitude of change can be calculated between groups of samples for a gene of interest.

However, this method of analysis is greatly complicated by the fact that housekeeping genes commonly used as normalizers (e.g., GAPDH, β -actin, and HPRT) can change expression levels when comparing tissues or cells in different states (Bustin, 2000; Schmittgen *et al.*, 2000; Goidin *et al.*, 2001; Hamalainen *et al.*, 2001). 18S rRNA is another normalizer that intuitively and experimentally appears to be more stable, but even 18S can vary in comparison to other genes when analyzed by sensitive real-time PCR techniques (Bustin, 2000, Akilesh *et al.*, 2003). Any small variation in the normalizer amplification would therefore compromise the analysis of the complete real-time PCR data set.

Global Pattern Recognition™ Analysis Tool processes the data derived from groups of biological replicates and reveals the changes in groups of genes. Global Pattern Recognition™ Analysis Tool positions the expression level of each gene with respect to all genes within an experiment. This can be done without prior assumption that a gene (normalizer) has an invariant expression level. Global Pattern Recognition™ Analysis Tool is unbiased in that it enables the experimental data, not the experimenter, to define the invariant normalizer genes. The use of any gene as a potential normalizer also maximizes the use

of the limited real estate on a StellarArray™ qPCR Array by eliminating the loss of wells used to contain potentially erroneously predefined normalizers.

3.0 Procedure Overview

Statistical data analysis requires the use of at least 3 biological replicates. Therefore, GPR uses a minimum of 3 'Control' samples versus 3 'Experimental' samples in order to facilitate a statistical analysis of the gene expression data. (Figure 1)

A typical experiment would utilize biological replicates. Biological replicates are defined as samples collected from separate and closely matched biological samples. They are processed individually under closely matched conditions. We recommend using a minimum of 3 biological replicates in each of two groups, for example '3 sick vs. 3 healthy' or '3 treated vs. 3 untreated' groups.





Figure 1. Example of experimental design to produce statistically correct data.

4.0 StellARray™ qPCR Protocol at a Glance



5.0 Plate Format

StellARray™ Plates are available in four different plate formats, each tailored to a specific subset of real-time PCR instruments and associated blocks.

Product type	For Real-Time Instruments
96-well Plate	Most thermal cyclers (non-ABI®, non-Roche®, non-low profile), such as Eppendorf®, Bio-Rad®, iCycler®, iQ5®, myiQ® Stratagene™, MX3000P® and MX3005P®
FAST 96-well Plate	ABI® FAST 96-well blocks
AB 96-well Plate	All ABI® "standard" blocks (7000, 7300, 7500, 7700, 7900)
384-well Plate	All common 384-well thermal cyclers, such as ABI®, Eppendorf®, Bio-Rad®

If you are using the Roche® LightCycler® 480 with 96-well block or 384-well block or a device requiring a low-profile plate such as the BIO-RAD® CFX96® or Stratagene™ MX4000®, please contact us to inquire about a StellARray™ qPCR Array in compatible plate format.

Note

Be sure that you have the correct PCR array format for your instrument before starting the experiment.

6.0 Materials Needed (not provided)

- Optical adhesive plate covers
- 2 x SYBR® Green Master Mix, with hot start *Taq* DNA polymerase

Note

Depending on the requirements of your thermal cycler, StellARray™ Products can be used with SYBR® Green Master Mixes containing ROX™, fluorescein or no passive reference dye. Check with your instrument manufacturer for the requirements regarding the passive reference dyes used in the Master Mix. Be sure to use the correct Master Mix for your thermal cycler before continuing with this protocol.

- Nuclease-Free dH₂O
- cDNA samples for gene expression analysis or gDNA samples for copy number variation (CNV) analysis.

Note

Ensure you have read the recommendations for sample preparation and quality control.

7.0 Recommendations for RNA Sample Preparation and Quality Control

Note

High quality RNA is ESSENTIAL for obtaining accurate gene expression results. The most important prerequisite for any gene expression analysis experiment is consistent, high-quality RNA from every experimental sample. Therefore, the sample handling and RNA isolation procedures are critical to the success of the experiment. Residual traces of proteins, salts and especially genomic DNA, will negatively impact the accuracy of gene expression experiments.

The use of high quality products for RNA purification is strongly recommended (please contact scientific.support@lonza.com for specific recommendations).

A DNase treatment step (on-column or in separate tube) must be performed.

RNA Quality Control

Remove an aliquot of the DNase-treated RNA for quantification and qualification. Quantification and check of purity should be done by UV spectrophotometry.

- $A_{260}:A_{230}$ ratio should be greater than 1.7

- $A_{260}:A_{280}$ ratio should be greater than 2.0
- Concentration by A_{260} should be greater than 4 $\mu\text{g}/\text{ml}$ total RNA

Determine the quality of your RNA samples by denaturing agarose gel electrophoresis (e.g. by using the FlashGel® System for RNA; Lonza)

Note

There should be a sharp distinction at the small side of both the 18S and 28S ribosomal RNA (rRNA) bands or peaks. Any smearing or shoulder to the rRNA bands or peaks indicates that degradation has occurred in the RNA sample.

8.0 cDNA Synthesis

A cDNA generation kit that works with random decamers (hexamers are acceptable) or oligo-dT primers or a mixture of both should be used for the cDNA synthesis step.

9.0 Protocol for qPCR

1. Calculate experimental cocktail volumes using the Reagent Calculator. (<http://array.lonza.com/reagentcalculator>)

Super Mix Recipe (per plate)

Note

At least 3 replicates are required per experimental condition

	96-well plate	96-well fast plate	384-well plate
2xSYBR® Green Master Mix	1056 μl	528 μl	2112 μl
dH ₂ O	998 μl	499 μl	1996 μl

Sample-Specific Working Solutions (per plate)

	96-well plate	96-well fast plate	384-well plate
Super Mix	2006 μ l	1003 μ l	4013 μ l
Template*	106 μ l	53 μ l	211 μ l

* Typically a 1:10 or 1:20 dilution of cDNA in nuclease-free water; gDNA 400 ng per 96-well plate; 1500 ng 384-well plate.

- Remove foil seal from StellarArray™ Plate.
- Either:** Add 20 μ l Sample-Specific Working Solution to each well (96-well plate)
Or: Add 10 μ l Sample-Specific Working Solution to each well for 384-well or 96-well fast plates.
- Firmly seal the StellarArray™ Plate.
- Centrifuge (hand-spin) the StellarArray™ Plate to ensure all reagents are at bottom of well.
- “Soak” StellarArray™ Plate by incubating at room temperature (or 4°C) for 15 minutes prior to loading into qPCR machine.
- Load plate into preprogrammed machine.

Standard PCR conditions

One (1) cycle at 50°C for 2 minutes

One (1) cycle of recommended Master Mix Specific Hot Start Temperature and Time

Forty (40) cycles at 95°C for 15 seconds followed by 60°C 1 for minute
Dissociation (Melt) Curve – refer to instrument specifications

Fast PCR conditions

One (1) cycle at 50°C for 1 minute

One (1) cycle of recommended Master Mix Specific Hot Start Temperature and Time

Forty (40) cycles at 95°C for 1 second followed by 60°C at 20 seconds
Dissociation (Melt) Curve – refer to instrument specifications

Perform Global Pattern Recognition™ Data Analysis

8. Prepare qPCR data for export to GPR. On threshold based instruments, set common threshold value (typically 0.15-0.2) after evaluating amplification curves.
9. Export data produced to a CSV (Comma Separated Value) or Microsoft® Excel® file spreadsheet. The format of the file should adhere to the following guidelines (Fig. 2):
 - The first field should contain well positions
 - The next group of fields should be the CT values from the control
 - The next group of fields should be the CT values from the experiment

WELL POSITION

	A	B	C	D	E	F	G	H	I	J	K
1	A01	34.17791	34.13614	34.42684	33.31817	31.61491	34.61817	33.53541	34.46147	36.62733	33.03477
2	A02	33.70426	31.77424	31.54926	30.47718	31.64141	31.90771	32.70826	32.08391	31.73788	31.42541
3	A03	29.33333	28.33333	28.33333	28.33333	28.33333	28.33333	28.33333	28.33333	28.33333	28.33333
4	A04	31.28889	31.28889	31.28889	31.28889	31.28889	31.28889	31.28889	31.28889	31.28889	31.28889
5	A05	32.88889	32.88889	32.88889	32.88889	32.88889	32.88889	32.88889	32.88889	32.88889	32.88889
6	A06	33.88889	33.88889	33.88889	33.88889	33.88889	33.88889	33.88889	33.88889	33.88889	33.88889
7	A07	34.88889	34.88889	34.88889	34.88889	34.88889	34.88889	34.88889	34.88889	34.88889	34.88889
8	A08	35.88889	35.88889	35.88889	35.88889	35.88889	35.88889	35.88889	35.88889	35.88889	35.88889
9	A09	36.88889	36.88889	36.88889	36.88889	36.88889	36.88889	36.88889	36.88889	36.88889	36.88889
10	A10	37.88889	37.88889	37.88889	37.88889	37.88889	37.88889	37.88889	37.88889	37.88889	37.88889

Figure 2. Example of layout for uploading data to Global Pattern Recognition™ Analysis Tool.

10. Upload data to Global Pattern Recognition™ Analysis Tool (<http://array.lonza.com/gpr>):
 - Enter the Lonza product identifier (For Example: 00188162)
 - Please select either a CSV (Comma Separated Value) or Microsoft® Excel® file to upload
 - Indicate the number of control samples and experiment samples

10.0 Interpretation of Results

The results obtained from the GPR analysis are not ranked strictly by fold change. When a single gene normalizer is selected, gene expression changes are typically ranked by their magnitude of change using the $\Delta\Delta Ct$ method, with those genes showing the largest fold changes ranked as most significant. Unfortunately, these large changes in gene expression may mask small but biologically important changes in gene expression, such as master regulator genes (e.g., transcription factors).

In biological systems, however, large is not always synonymous with importance. The Global Pattern Recognition™ Analysis Tool algorithm generates a ranked list of significantly changed genes within a real-time PCR dataset. The unique algorithm and accompanying software overcomes the problem of identifying invariant normalizers and the pitfalls of producing faulty statistics based merely on magnitude of change. Only after the genes are statistically ranked is the magnitude of the change calculated. The Global

Rank	Gene Name	p-value	Fold Change	Rank	Gene Name	p-value	Fold Change
1	FOS	0.0001	29.7582	20	SLC16A3	0.0240	1.6755
2	EGR1	0.0004	20.6085	21	FADD	0.0242	-1.8505
3	JUNB	0.0005	7.9065	22	ABCC2	0.0259	-1.5363
4	GADD45B	0.0015	4.0261	23	SLC22A6	0.0263	5.1732
5	CYP1A1	0.0023	-3.5846	24	TXNRD1	0.0270	-1.5292
6	EDN1	0.0046	-2.4949	25	SLCO1C1	0.0276	3.6201
7	HMOX1	0.0049	-2.6288	26	MMP9	0.0280	2.6985
8	MYC	0.0051	-2.6534	27	TNF	0.0308	1.9693
9	JUN	0.0056	2.3803	28	NAT2	0.0347	-1.5462
10	DNAJB1	0.0062	2.3285	29	CCNG2	0.0378	1.6557
11	HBEGF	0.0067	3.7693	30	ABCB9	0.0381	-1.6587
12	SERPINE1	0.0081	3.0827	31	VEGFA	0.0390	1.6334
13	CYP3A4	0.0091	-2.4592	32	BCL6	0.0396	1.6786
14	IGFBP3	0.0131	2.0451	33	ESR1	0.0398	1.8955
15	NOS2A	0.0137	2.7647	34	BCL2L1	0.0400	-1.4412
16	MCL1	0.0140	-1.8457	35	FOSB	0.0404	2.4147
17	GAP43	0.0153	3.1655	36	CYP7A1	0.0405	2.6760
18	MT3	0.0160	2.5987	37	SLC2A1	0.0427	1.6791
19	CYR61	0.0177	1.8896				

Figure 3. GPR results derived from Acetaminophen treatment of Clonetics® ReadyHeps™ Primary Human Hepatocytes. Results shown represent genes whose expression levels are statistically different between the treated and untreated groups and Fold Change values are displayed with respect to the treated group.

Pattern Recognition™ Tool provides a true statistical analysis of results based on consistency in the data, making the Global Pattern Recognition™ Analysis Tool optimally suited to detect small, but reproducible changes.

11.0 References

1. Wo Y, Wright SM, Maas SA, Alley TL, Caddle LB, Kadmar S, Affourti J, Foreman O, Shaffer D, Bronson RT, Roopenian DC, Mills KD. The nonhomologous end joining factor Artemis suppresses multi-tissue tumor formation and prevents loss of heterozygosity. *Oncogene*. 2007 Sep 6;26(41):6010-20.
2. Seymour RE, Hasham MG, Cox GA, Shultz LD, Hogenesch H, Roopenian DC, Sundberg JP. Spontaneous mutations in the mouse Sharpin gene result in multiorgan inflammation, immune system dysregulation and dermatitis. *Genes Immun*. 2007 Jul;8(5):416-21.
3. Ueno M, Lyons BL, Burzenski LM, Gott B, Shaffer DJ, Roopenian DC, Shultz LD. Accelerated wound healing of alkali-burned corneas in MRL mice is associated with a reduced inflammatory signature. *Invest Ophthalmol Vis Sci*. 2005 Nov;46(11):4097-106.
4. Brown AC, Olver WI, Donnelly CJ, May ME, Naggert JK, Shaffer DJ, Roopenian DC. Searching QTL by gene expression: analysis of diabetes. *BMC Genet*. 2005 Mar 10;6(1):12.
5. Chen YG, Choisy-Rossi CM, Holl TM, Chapman HD, Besra GS, Porcelli SA, Shaffer DJ, Roopenian D, Wilson SB, Serreze DV. Activated NKT cells inhibit autoimmune diabetes through tolerogenic recruitment of dendritic cells to pancreatic lymph nodes. *J Immunol*. 2005 Feb 1;174(3):1196-204.

12.0 FAQs

What type of Master Mix should we use with the StellARray™ Gene Expression System?

The Master Mix of choice is a 2x SYBR® Green Master Mix with a hotstart *taq* polymerase. Check with the instrument manufacturer for the requirements regarding the passive reference dyes used in the Master Mix. StellARray™ qPCR Arrays are compatible with Master Mixes containing ROX™ (e.g. for ABI® and Stratagene® Thermal Cyclers), fluorescein (e.g. for Bio-Rad® iCycler®, MyiO™, and iQ™5 Thermal Cyclers) or without (e.g. Roche®, Bio-Rad® Opticon™, Opticon™ 2, and Chromo4™ Thermal Cyclers) passive reference dye.

How are genes ranked in the GeneSieve™ Query?

All mammalian genes tagged with the search term are filtered through the GeneSieve™ Query modules to look for relationships. The result is a hit list of genes ranked by their relevance associated with the particular search term. GeneSieve™ Query then filters genes of interest using a variety of parameters such as overall score, the type of data used, literature reference, and microarray fold change.

For gene expression, are the primers designed to detect 3'-end or 5'-end mRNA sequences?

The primer designs are not usually directed to any specific region of a message as the aim is for primer candidates that will produce highly efficient and specific PCRs. Because the use of a mixture of oligo-DT and random decamers is recommended, the reverse transcription is not effected by the location of the amplicon.

If the same assay can detect both gDNA and mRNA, how do I know that your plate won't detect genomic DNA when I am trying to measure changes in mRNA expression?

As for any method for gene expression analysis, quality of RNA sample prep is critical. StellarRay™ Products can be used to detect genomic DNA or cDNA (synthesized from mRNA). If you are evaluating mRNA expression, genomic DNA must be removed by DNase digestion before cDNA synthesis. StellarRay™ Plates include a control that checks for the absence of genomic DNA.

Does the StellarRay™ Process and GPR analysis use normalizers?

Yes, but not in the "conventional" sense. GPR finds the genes that are invariant and sets those as "normalizers". The researcher no longer has to guess which housekeeping gene is invariant before conducting the experiment

How much RNA is needed per well in a 96-well plate?

250 to 500 ng total RNA is recommended per cDNA reaction. A dilution of the cDNA is made (typically 1:20 in dH₂O) and the required volume of cDNA

added to the Sample-Specific Working Solutions. In some cases successful experiments have been conducted using as little as 25 ng total RNA per cDNA reaction and applied to a 384-well StellarArray™ qPCR Array.

How are the primers validated?

Each primer pair undergoes extensive *in silico* and wet lab validations, ensuring specificity and high sensitivity. Primer design parameters were selected to enhance:

- Gene specificity
- Amplicon efficiency
- Uniformity of primer Tm's
- A/T-rich 3'-ends of primers
- Elimination of primer-dimer artifacts
- Generation of single-peak dissociation curves
- Detection of single bands of expected size via gel electrophoresis

By utilizing stringent criteria for both design and performance, accurate, gene-specific results are ensured.

If technical support is required please contact
scientific.support@lonza.com

13.0 Ordering Information

Please find the ideal StellarArray™ qPCR Array for your research area at www.lonza.com/arrays.

14.0 Related Products

FlashGel® System for RNA

Contact Information

North America

Customer Service: 800 638 8174 (toll free)
Scientific Support: 800 521 0390 (toll-free)
scientific.support@lonza.com
Online Ordering: www.lonza.com

Europe

Customer Service: +32 87 321 611
Cell Discovery
Scientific Support: +49 221 99199 400
scientific.support.eu@lonza.com
Molecular Biology, RTS & Media
Scientific Support: +32 87 321 611
scientific.support.eu@lonza.com
Online Ordering: www.lonza.com

International

Contact your local Lonza Distributor
Customer Service: +1 301 898 7025, ext. 2322
scientific.support@lonza.com

International Offices

Australia +61 3 9550 0883
Austria 0800 201 538 (toll free)
Belgium +32 87 321 611
Brazil +55 11 2069 8800
Denmark +45 43 56 74 00
France 0800 91 19 81 (toll free)
Germany 0800 182 52 87 (toll free)
India +91 22 4342 4000
Italy +39 0363 45710
Japan +81 3 5566 0612
Poland +48 22 833 87 45
Singapore +65 64914214
Spain +34 902 531 366
Sweden 020 140 4410 (toll free)
Switzerland 0800 83 86 20 (toll free)
The Netherlands 0800 022 4525 (toll free)
United Kingdom +44 118 979 5234

Lonza Walkersville, Inc.
Walkersville, MD 21793

For Research Use Only. Not for use in diagnostic procedures.

GeneSieve, StellARray, and Global Pattern Recognition are trademarks of Bar Harbor BioTechnology, Inc., and are manufactured exclusively for Lonza. Unless otherwise noted, all trademarks herein are marks of the Lonza Group or its affiliates. ABI is a registered trademark of Applied Biosystems Inc.; Roche and LightCycler are registered trademarks of Roche Inc.; Eppendorf is a registered trademark of Eppendorf-Netheler-Hinz GmbH; Stratagene is a trademark of Agilent Technologies; SYBR is a registered trademark of Molecular Probes Inc.; Bio-Rad, iCycler, Mj10, iQ5 Opticon and Chromo4 are trademarks of Bio-Rad Laboratories, Inc.; ROX is a trademark of Applied Corporation; Microsoft and Excel are registered trademarks of Microsoft Corporation

The information contained herein is believed to be correct and corresponds to the latest state of scientific and technical knowledge. However, no warranty is made, either expressed or implied, regarding its accuracy or the results to be obtained from the use of such information and no warranty is expressed or implied concerning the use of these products. The buyer assumes all risks of use and/or handling. No statement is intended or should be construed as a recommendation to infringe any existing patent.

© Copyright 2009, Lonza Walkersville, Inc.
All rights reserved.
MN-StellARrayGeneExpress 04/09
